



MRX001



M E R X I N TM

MRX001 is a generic blister multidose dry powder inhaler

**It is an off-the-shelf pre-assembled device
ready to be filled with your formulations**

MRX001

MRX001 is a generic blister based multidose dry powder inhaler (DPI) device.

MRX001 is based on the open-click-inhale-close principle and can be used as a generic device of current branded DPIs.

It is supplied as a pre-assembled device ready to be filled with a blister strip manufactured at client's facilities in the country of their choice, to meet local requirements.

MRX001 uses Merxin's proprietary designs and does not knowingly infringe any IP.

MRX001 is available with all the necessary documentation for regulated markets such as EU, US and Japan as well as emerging markets.

Why work with Merxin?:


You will save time and money in using a readymade first class device.

You will access a popular and established device technology.

You will tap into 1st class expertise on quality documentation ready for EU and US markets.

MRX001 features:

- A multidose dry powder inhaler blister device (blister not supplied in standard configuration)
- Up to 60 doses
- Dose counter
- Substitutable device
- Open-click-inhale-close
- Overall dimensions: Diameter 80 mm x thickness 27 mm
- Pharmaceutical grade materials for inhalation devices (ABS, POM, PP)
- Customisable colours
- Proprietary design



About Merxin:

Merxin is a generic medical device development company.

Development and supply of generic inhalation devices as "off the shelf and ready to use" for finished product development is our core business.

Merxin's portfolio is based on a generic device of a multidose blister dry powder inhaler, open-click-inhale-close.

By supplying standardised off-the-shelf devices, Merxin provides clients with a lower cost and lower risk route to market.

Combined with the expertise and support of the Merxin team, you can develop and deliver your products easily and efficiently.

Get in touch to ask about terms of access and prices



M E R X I N™

Contact:

Merxin can be contacted at info@merxin.com and via our website: www.merxin.com

Generic Medical Devices Made Easy