


**MRX002**



M E R X I N <sup>TM</sup>

MRX002 is a clear PET Vial to study HFA formulations

## Be in touch now to expand and develop your inhaled portfolio

### MRX002:

MRX002 is a PET vial ideal for modelling pMDI formulations.

MRX002 is a clear PET vial that can be fitted with commercial pMDI valves, continuous and metered.

With MRX002, formulations can be observed and understood in situ. MRX002 is safe to use with fluorinated liquids and a range of solvents compatible with PET. Note however, that MRX002 is not yet pressure rated and is for laboratory use only.

MRX002 uses Merxin's proprietary designs and does not knowingly infringe any IP.

In addition to supplying MRX002, Merxin can support clients with formulating, testing, manufacturing, and clinical and registration expertise via its partner organisations:

Aedestra  
[www.aedestra.com](http://www.aedestra.com)

and

Lustren  
[www.lustren.com](http://www.lustren.com)

### MRX002 features:

-Clear PET vial to study model pMDI formulations in situ.

-Dimensions:  
Height 75 mm  
Diameter 22 mm  
Crimp diameter 20 mm

-Material: PET (Polyethylene terephthalate).

-Compatible with HFA solvents, ethanol, water, PEG and glycerol.

-Example of compatible pMDI valves:

VARI KHFA range  
Bespak MDI valves

-Available in quantities of your choice.

-Ideal for formulation screening. For laboratory use only.



### About Merxin:

Merxin is a generic medical device development company.

Development and supply of generic inhalation devices as "off the shelf and ready to use" for finished product development is our core business.

By supplying standardised off-the-shelf devices, Merxin provides clients with a lower cost and lower risk route to market.

Combined with the expertise and support of the Merxin team, you can develop and deliver your products easily and efficiently.

**Get in touch to ask about terms of access and prices**



Generic Medical Devices Made Easy

### Contact:

Merxin can be contacted at [info@merxin.com](mailto:info@merxin.com) and via our website: [www.merxin.com](http://www.merxin.com)